

IMPROVING BOND STRENGTH THROUGH ACID ETCHING OF DENTIN AND BONDING TO WET DENTIN SURFACES

JOHN KANCA III, D.M.D.

Restorative dentistry needs a dentin-enamel bonding system that can:

- provide gapfree composite restorations;
- achieve rapidly developing high-bond strengths;
- be biocompatible;
- function in the presence of moisture;
- treat dentin and enamel simultaneously;
- treat multiple types of intraoral surfaces.

Producing a gapfree restoration is critical, especially in composite resin restorations. Brannstrom has shown that bacteria can rapidly colonize the gaps between composite resin restorations and the tooth interface.¹ Bergenholtz and others found that bacteria are directly correlated to pulpal pathosis.²

Bonding systems need to generate high-bond strengths in a short time. Typically, 24-hour bond strengths are the reported values, but no one waits 24 hours to begin inserting composite resins. Davidson and others showed that the polymerization shrinkage of composite resins places stress on bonding systems.³ To resist the contraction force of polymerizing composites, bond strengths must develop rapidly. The need for biocompatibility is

ABSTRACT

A bonding system using moisture on the tooth surface can be an enormous benefit as obtaining dentin dryness in the mouth is nearly impossible. The author describes a system that bonds to both wet and dry surfaces.

obvious, but it's likely that all current bonding systems are biocompatible.

Since dentin is inherently a wet substrate, the bonding system must function in the presence of moisture. The need for a common treatment of dentin and enamel is again obvious to the practicing dentist, who well knows there really is no other way to treat. Often it's very difficult to recognize where enamel ends and dentin begins. Keeping an enamel etchant off dentin all the time is not practical. Remember this when considering a bonding system that is detrimentally affected by the enamel etchant contacting the dentin surface.⁴

Finally it would be most convenient to the practitioner if one bonding system could be

used to bond all types of surfaces. Rather than purchase one kit to bond enamel, another for metal and yet another for porcelain, it would be highly desirable to be able to use the same bonding system for all surfaces with only small alterations. For this reason we should consider not dentin enamel bonding systems but rather adhesive systems.

Perhaps the most intriguing need is function on moist or wet dentin surfaces. Freshly cut dentin presents a difficult restorative medium for adhesive resin systems. The dentin substrate is inherently moist because of the tubular fluid. Reports have demonstrated that moisture typically deters the ability of resin bonding systems to adhere to dentin. Terkla and others⁵ showed that two dentin bonding resin systems tested could not seal moist dentin.

Mitchem and others⁶ noted that the bond strengths of Scotchbond (3M Dental Products), Tenure (Den-Mat) and Gluma (Columbus Dental) fell when tested under simulated physiological conditions. Andreaus and others⁷ reported that water in dentinal tubules interfered with the adhesion of dentin adhesives. Tao and Pashley⁸ reported that the bond



Figure 1. Maxillary left central incisor before restoration of dentin.



Figure 2. Cord is packed in gingival sulcus to prevent fluid contamination.

strength of Gluma and Scotchbond fell significantly when smear layers were removed and pulpal pressure was applied to a dentin surface. Prati and others⁹ found that storage in simulated pulpal pressure caused a significant reduction in the bond strength of Clearfil Photobond (J. Morita) and Scotchbond 2 in deep dentin. Mitchem and Gronas¹⁰ and Glasspoole and others¹¹ reported that the bond strength of Scotchbond 2 was diminished in the presence of dentinal

moisture.

This report presents data demonstrating the improvement in bond strengths obtained with a particular bonding system when dentin is acid-etched and the primers are applied to wet dentin surfaces. The clinical technique is also illustrated.

MATERIALS AND METHODS

Sixty extracted human molars were selected for this study. They were free of obvious

defects and used within two weeks of harvest. They were cleaned with a soap containing chlorhexidine and placed in a 3 percent glutaraldehyde solution for 24 hours. They were then rinsed and stored in tap water at 4 C until time of usage. For the experiment the teeth were embedded in an autopolymerizing acrylic resin in a stainless ring with an inside diameter of one inch. The teeth were ground on a model trimmer to expose adequate facial dentin surface for bonding.

The samples were stored in water at 37 C for 24 hours. After storage, the exposed dentin surfaces were wet-sanded with 320 grit sandpaper just before treatment. The teeth were randomly divided into three equal groups, each of which was divided into two subgroups.

The dentin-enamel bonding system examined was the All Bond 2 system (Bisco Dental). The system has available two concentrations of phosphoric acid semigel etchant conditioners: a 10 percent and a 32 percent concentration. The system also contains a Primer A and Primer B and a dual-cure capability unfilled bonding resin. Primer A is made up of NTG-GMA(N-tolyl glycine

TABLE 1

LIST OF MATERIALS

All-Etch (10% phosphoric acid semigel etchant)
Uni-Etch (32% phosphoric acid semigel etchant)
Primer A (batch 129161)
Primer B (batch 129171)
Dentin-enamel bond (batch 129091)

glyc-
and
(bipl
acet
T.
thre
trea
man
— G
was
surfi
drie
of th
cent
surfi
A an
toge
mixt
dent
five
surfi
coat
E.
spre
the
fifth
surfi
for 5
drie
the
bond
appl
air-t
and
seco
light
Rese
milli
cylir
dian
the
mat
dent
with
purp
shad
comj
for 4
O.
sme
visib
mois
in th

glycidylmethacrylate)/acetone and Primer B is BPDM (biphenyl dimethacrylate)/acetone (Table 1).

The teeth were divided into three equal groups and were treated in the following manner:

■ Group 1. The smear layer was not removed. Dentin surfaces in 10 teeth were air-dried for 3 seconds with the tip of the air syringe about 3 centimeters from the tooth surface. A drop each of Primer A and Primer B were mixed together in a well. The primer mixture was applied to the dentin surface in a minimum of five consecutive coats. The surfaces were not dried between coats.

Each application was able to spread out on the surface before the next was made. After the fifth application, the primed surface remained undisturbed for 5 seconds and then was air-dried for 5 seconds. A layer of the dentin-enamel unfilled bonding resin (Bisco) was applied to the primed surface, air-thinned for 1 to 2 seconds and light-polymerized for 20 seconds with a Demetron 401 light-curing unit (Demetron Research). A Teflon mold 2.5 millimeters thick with a cylindrical matrix 4 mm in diameter was then clamped to the embedded tooth so that the matrix was over the treated dentin. The matrix was filled with a single increment of All-purpose composite (universal shade, Bisco Dental), and the composite was light-activated for 40 seconds.

On the second 10 teeth, the smear layer was rinsed but left visibly wet, or moist. The moisture presence was obvious, in that the dentin surface was




Figure 3. Cavity is prepared with light bur treatment to enhance retention.




Figure 4. Preparation is treated with semigel etchant and rinsed with air-water.

shiny with moisture. A facial tissue was sprayed with air-water until it was wet and wiped the surface of the samples, removing only the excess water. (For a patient, simply use a disposable brush to rewet the dentin surface should it become dry.) A drop each of Primer A and B were mixed together in a well and applied to the wet dentin surface in a minimum of five consecutive coats as described. The primed surface was dried

and the dentin-enamel bond applied. Composite cylinders were bonded to the treated dentin surface.

■ Group 2. The dentin was treated with 10 percent phosphoric acid gel etchant (All-Etch, Bisco Dental Products) for 30 seconds and rinsed for 5 seconds with an air-water mixture. On 10 teeth, the dentin was dried for 3 seconds with an air syringe at full blast at a distance of about 2 cm. On the remaining 10 teeth, the



Figure 5. Matte appearance of surface shows area too dry to be treated.

dentin was left wet as described in group 1. Primer A and B and the dentin-enamel bonding resin were applied as described. It sometimes took more than five coats of the Primer A and B mixture to produce a glossy coating of dentin. We placed enough applications of the A and B mixture to achieve this appearance. Composite cylinders were bonded to the dentin surfaces as previously described.

Group 3. The dentin in this group was treated with 32 percent phosphoric acid gel etchant (Uni-Etch, Bisco Dental) for 20 seconds and rinsed for 5 seconds with an air-

water mixture. On 10 teeth, the dentin was dried for 3 seconds with an air syringe as in group 2. On the remaining 10, we wiped the dentin with a wet tissue to remove only the excess water as in group 1. Primer A and B and the dentin-enamel unfilled bonding resin were applied to the dentin as described in group 2. Composite cylinders were bonded to the samples as above.

For a control group, 10 incisors were embedded in acrylic and ground to produce adequate enamel surface for bonding. The enamel was treated with the 32 percent etchant for 20 seconds, rinsed

and dried. Primer A and B and the bonding resin were applied as in group 3. Composite cylinders were applied as previously described.

The matrices were removed from the teeth and the samples were stored for 24 hours in water at 37 C. Shear bond strength was measured on an Instron machine (Model no. 1123) with a crosshead speed of 5 mm/minute and recorded. We used a sharp blade to shear the samples with the blade parallel and immediately adjacent to the bonded dentin surface.

RESULTS

The best results were obtained when the dentin was acid-etched and the primer mixture applied to wet dentin surfaces (Table 2). The lowest bond strengths occurred when the primer mixture was applied to a dry smear layer. In all cases, we obtained higher bond strengths when the primer mixture was applied to a wet dentin surface. Failure modes were examined with a light microscope.

In the 10 percent etch/wet and 32 percent etch/wet groups, all samples failed cohesively within the body of the dentin. In the 10 percent etch/dry and 32 percent etch/dry groups, four of 10 samples failed cohesively in the dentin. All samples bonded to a smear layer experienced adhesive failure between the resin and the smear layer. Two-way ANOVA analysis of each of the

TABLE 2

BOND STRENGTHS WITH VARIOUS SURFACE TREATMENTS

GROUP	SURFACE DRY			WET	
	Mpa (S.D.)	CV	Mpa (S.D.)	CV	
320 grit smear layer	11.7 (4.8)	46.3%	24.1 (2.1)	8.8%	
10% H3PO4 30s	24.3 (5.2)	21.4%	34.3 (2.4)	7.0%	
32% H3PO4 20s	20.7 (10.8)	52.2%	36.5 (1.6)	4.4%	
32% H3PO4 20s (enamel)	31.9 (4.8)	15.0%			

All groups N = 10.
CV = coefficient of variance.

Figur

Figur

grou
dry)
bonc
sign
thar
(P<
vari
for
sam
gree

Dis
The
vita
tubi
moi



Figure 6. Preparation surface must be wet for higher bond strength.




Figure 7. Primers are applied directly on wet dentin.

groups (etch vs. no etch, wet vs. dry) shows that in each case the bond strength of materials was significantly higher to the wet than the dry dentin surface ($P < .05$). The coefficient of variation, however, was lower for all of the wet-bonded samples. This may indicate a greater reliability of the bond.

DISCUSSION

The natural state of dentin in vital teeth includes fluid in the tubules. This tends to create a moist surface on dentin.

Nordenvall reported that it is difficult to dry vital dentin and to maintain a dry state because of the outward flow of dentinal fluid.¹² A dentin-enamel bonding system that can adhere to a wet surface has a great advantage. It's difficult to determine that the tooth surface is completely dry before applying conventional dentin bonding materials that demand dryness. This is particularly true in posterior areas and any areas to which access is compromised or difficult.

The most intriguing aspect of the material being studied is that we obtained the best adhesion values with wet surfaces. As we explained in an earlier report, mixing the acetone in the primer mixture with the water on the dentin surface caused the vapor pressure of the water to rise and created some evaporation.¹³ Mixing the acetone/resin mixture with the water also reduces the surface tension of the water so that the mixture seems to "chase" the water.

This results in the primer layer being laid down in intimate contact with the dentin surface; that is, the intertubular dentin, peritubular dentin and presumably the dentinal tubule lumina are adapted to more thoroughly when the primer mixture is applied to the wet surface. This system can generate gapfree restorations in vitro and in vivo, and the resin penetration in vivo was the same as in the in vitro state.¹⁴ (In this report, wetness or moisture refers only to dentinal fluid, water or a disinfectant solution. Salivary contamination is unacceptable as saliva contains viscous proteinaceous material as well as bacteria.)

It has also been reported that the first version of the All Bond system can bond to wet enamel as well as to wet dentin.¹⁵ Bond strengths obtained in that report were on the order of 29 to 32 Mpa. The current data show this system's ability to generate bond strengths to dentin in excess of those to enamel. (Data to be published show this to be the same with the All Bond 2 system.)

All discussions of dentin bonding materials should contain some reference to an




Figure 8. Five coats of primer are applied without drying. Preparation is gently air-dried.



Figure 9. The dentin-enamel bonding resin is applied to primed surface.

etched enamel standard. The simple disclosure of dentin bond strength values without that standard can be difficult to interpret. This bonding system uses a phosphoric acid dentin-enamel conditioner. Typical reports on dentin bonding cite two references that noted that applying phosphoric acid to dentinal surfaces reduced bond strengths.^{16,17} The authors used a hydrophobic unfilled BIS-GMA autopolymerizing resin on etched dentin surfaces with 50

percent phosphoric acid on the dentin for up to 2 minutes.

Both the concentration of acid and time of application have changed considerably. It has still been reported, however, that the bond strengths of even current bonding systems may be diminished by using phosphoric acid as the dentin conditioner.⁴ The greatest controversy in this system stems from the belief that applying phosphoric acid to cut dentin will harm the pulp. This

was concluded in some early studies examining the alleged effects of acid on dentinal surfaces.^{18,19}

In those studies, Class V cavity preparations were cut into teeth and the cavities in the control group were restored with zinc oxide-eugenol. In the experimental group, the cavities were acid-etched with 50 percent phosphoric acid for 1 minute and then rinsed, dried and filled with ZOE. The pulps in the experimental groups showed a negative reaction to something, and the authors concluded that that something was the acid treatment.

But a review of the literature showed that the control material (ZOE) is cytotoxic when placed against living tissue.²⁰⁻²⁷ In addition, acid treatment of dentin increases its permeability.²⁸ Eugenol can diffuse through intact dentin,²⁹ and Brannstrom and coworkers have shown that ZOE-type materials have damaging effects on the pulp when placed on thin remaining dentin *in vivo*.²⁹⁻³² In an earlier report,³³ I speculated that the acid treatment of the dentin alone did not cause the pulpal reactions seen in those studies. It was more likely that ZOE application to the etched permeable dentin caused the reactions as eugenol (a phenol) diffused into the pulp more readily than when a smear layer was present, as in the unetched group. Considering these findings, we see that the major concern is not the acid treatment of dentin but the material placed on the dentin after etching and its ability to seal the dentin.

Recent reports support this contention. At the most recent AADR meeting, a paper showed

the A
biocc
after
treat
for u
react
less t
unet
prese
resul
with
and t
with
a ligh
findi
the fi
tive c

TEC

We r
centr
was t
the c
and t
shiny
faces
probl
than
recon
denti
bur to
using
produ
restor
Ste



Dr. Kanc
assistan
Universi
Alabama
ham and
Western
Universi
maintain
practice
restorati
390 Mid
Middleb
06782-2
Address
requests
author.

the All Bond bonding system as biocompatible when applied after 37 percent phosphoric acid treatment of freshly cut dentin for up to 2 minutes. The pulp reaction to this treatment was less than that to the ZOE unetched controls.³⁴ Another presentation noted the same results.³⁵ Dentin was treated with 37 percent phosphoric acid and the teeth were restored with New Bond (J. Morita) and a light-cured composite. These findings can significantly affect the future of adhesive restorative dentistry.


TECHNIQUE

We restored a maxillary left central incisor (Figure 1). There was tooth tissue absent from the cervical areas of the tooth, and the dentinal surfaces were shiny. These shiny dentin surfaces pose a difficult restorative problem as it's harder to restore than normal dentin.³⁶ We recommend preparing sclerotic dentin lightly with a diamond bur to enhance retention and using a bonding system to produce gapfree composite resin restorations in vivo.¹⁴

Steps for restoring this

maxillary left central incisor follow:

- Pack cord in the gingival sulcus to prevent fluids from contaminating the field (Figure 2).
- Prepare the cavity conservatively by light treatment with a diamond bur (Figure 3). No mechanical



Dr. Kanca is adjunct assistant professor, University of Alabama/Birmingham and Case Western Reserve University, and maintains a private practice in restorative dentistry, 390 Middlebury Rd., Middlebury, Conn. 06762-2634. Address reprint requests to the author.




Figure 10. Appearance after resin is applied and air-thinned.



Figure 11. Composite resin is applied incrementally and light-polymerized.

retention is used.

- Treat the preparation with the 10 percent semigel etchant for 30 seconds on both dentin and enamel and rinse thoroughly with air-water (Figure 4).
- If desired, use the 32 percent gel etchant on the dentin and enamel for 20 seconds in lieu of the 10 percent etchant. Figures 5 and 6 show the contrast between wet and dry prepared surfaces. As this bonding system obtains higher bond strengths to wet dentin, be sure

the preparation surface appears wet as shown in Figure 6. If the preparation surface becomes too dry, it will have a matte appearance (Figure 5). If this happens, remoisten the surface of the preparation with water or an antimicrobial solution (for example, Tubulicid Red, Global Dental Products) with a disposable brush.

- Mix together a drop of Primer A and Primer B in the well and place directly on the wet dentin (Figure 7).
- Apply five consecutive coats




Figure 12. Completed restoration at one-year recall.

(or more—use up the mixture) without drying between coats. After a few seconds, gently air-dry the surface for 2 to 3 seconds to remove any remaining volatiles (Figure 8).

■ Apply the dentin-enamel bonding resin to the primed surface and air-thin (Figure 9,10). Light-polymerize for 20 seconds.

■ Place composite resin (Herculite XR-V, Kerr Dental) incrementally and light-polymerize appropriately (Figure 11). Finish the filled preparation with burs and disks.

The completed restoration is shown in Figure 12. Adjacent to the newly completed restoration is a similarly restored tooth at one-year recall. I've used this system successfully in the past three years, retaining 58 of 58 Class V restorations over a 30-month period with no post-operative sensitivity. The All Bond system has also been shown to bond to all typical intraoral substrates.^{37,38}

CONCLUSION

This report shows that a dentin-enamel bonding system can not only bond to wet dentin sur-

faces, but also offer the best results when the dentin surface is wet. In addition, the best shear bond strength data are obtained through acid-etching the dentin. This provides a significant reduction in technique sensitivity. The clinician may apply the gel etchant to both the dentin and enamel at the same time—this is far easier than trying to etch the enamel alone. After rinsing the etchant from the tooth surface, remove only the excess water with a short blast of air—desiccation is not required.

High in vitro bond strengths result from this type of treatment, and the bonding system is biocompatible. It must be made clear that this report does not indicate endorsement of dentin treatment with any phosphoric acid and the indiscriminate placement of any bonding system. It has been reported that phosphoric acid treatment can diminish the bond strength of many bonding systems,⁴ and the manufacturers' instructions must be followed for each system. However, there are obvious advantages to using a system

that can treat dentin and enamel at the same time with phosphoric acid while improving the bond strength to both the enamel and the dentin. ■

Information about the products mentioned in this article may be available from the author. Neither the author nor the American Dental Association has any commercial interests in the products mentioned in this report.

1. Brannstrom M. Infection beneath composite resin restorations: can it be avoided? *Oper Dent* 1987;12:158-63.
2. Bergenholtz G, Cox C, Loesche W, Syed S. Bacterial leakage around dental restorations: its effect on the dental pulp. *J Oral Pathol* 1982;11:439-50.
3. Davidson C, DeGee A, Feilzer A. The competition between the composite-dentin bond strength and the polymerization contraction stress. *J Dent Res* 1984;63:1394.
4. Nathanson D, Amin F, Ashayeri N. Dentin etching vs. priming: effect on bond strength in-vitro. *J Dent Res* 1992;71(Special Issue):255(Abstract no. 1192).
5. Terkla L, Brown A, Hainisch A, Mitchem J. Testing sealing properties of restorative materials against moist dentin. *J Dent Res* 1987;66:1758-64.
6. Mitchem J, Terkla L, Gronas D. Bonding of resin dentin adhesives under simulated physiologic conditions. *Dent Mater* 1988;4:351-3.
7. Andreus S, Bayne S, Heymann H, Kanoy B. Intrapulpal composition fluid flow effects on dentin bond strength. *J Dent Res* 1989;(Special Issue):321(Abstract No 1114).
8. Tao L, Pashley D. Dentin perfusion effects on the shear bond strengths of bonding agents to dentin. *Dent Mater* 1989;5:181-4.
9. Prati C, Pashley D, Montanari G. Hydrostatic intrapulpal pressure and bond strength of dentin bonding systems. *Dent Mater* 1991;7:54-8.
10. Mitchem J, Gronas D. Adhesion to dentin with and without smear layer under varying degrees of wetness. *J Prosthet Dent* 1991;66:619-22.
11. Glasspoole E, Erickson R, Pashley D. The effect of surface moisture on bonding to dentin. *J Dent Res* 1991;70(Special Issue):457 (Abstract No. 1532).
12. Nordenvall K-J. Enamel bond-composite fillings in etched cavities: the appearance of resin surfaces facing enamel and dentine. *Swed Dent J* 1978;1:141-50.
13. Kanca J. Resin bonding to wet substrate. I. Dentin. *Quintessence Int* 1992;23:39-41.
14. Gwinnett AJ, Kanca J. Interfacial morphology of composite resin and shiny erosion lesions. *J Dent Res* 1992;71(Special Issue):133 (Abstract No. 224).
15. Kanca J. Resin bonding to wet substrate. II. Enamel. *Quintessence Int* 1992. In press.
16. Rider M, Tanner A, Kenny B. Investigation of adhesive properties of dental composite materials using an improved tensile test procedure and scanning electron microscopy. *J Dent Res* 1977;56:368-78.
17. Torney D. The retentive ability of acid-etched dentin. *J Prosthet Dent* 1978;39:169-72.
18. Macko D, Rutberg M, Langeland K. Pulpal response to the application of phosphoric acid to dentin. *Oral Surg Oral Med Oral Pathol* 1978;45:930-46.
19. Retief D, Austin J, Fatti L. Pulpal response to phosphoric acid. *J Oral Pathol* 1974;3:114-22.
20. Hensten-Pettersen A, Helgeland K. Evaluation of biologic effects of dental materials using four different cell culture techniques. *Scand J Dent Res* 1977;85:291-6.
21. Dahl B, Orstavik J. Cytotoxicity of temporary crown and bridge materials. *J Oral Rehabil* 1976;3:341-7.
22. Glass R, Zander H. Pulp healing. *J Dent Res* 1949;28:97-107.
23. Sayegh F, Reed A. Tissue reactions to a new restorative material. *J Prosthet Dent* 1969;22:468-78.
24. Spangberg L, Langeland K. Biologic effect of dental materials I. Toxicity of root canal filling materials on HeLa cells in vitro. *Oral Surg Oral Med*

Oral Pathol 1973;35:402-14.

25. Pascon E, Langeland K. Cytotoxicity of a new endodontic sealer. J Dent Res 1989;68(Special Issue):244(Abstract No. 501).

26. Kozam G, Mantell G. The effect of eugenol on oral mucous membranes. J Dent Res 1978;57:954-7.

27. Cook D, Taylor P. Tissue reactions to improved zinc oxide-cements. J Dent Child 1973;70(3):199-207.

28. Pashley D, Michelich V, Kehi T. Dentin permeability: effects of smear layer removal. J Prosthet Dent 1981;46:551-7.

29. Hume W. An analysis of the release and the diffusion through dentin of eugenol from zinc oxide-eugenol mixtures. J Dent Res 1984;63:881-4.

30. Brannstrom M, Billberg B. Pulp changes beneath temporary fillings with Pharmatec and zinc oxide-eugenol. Odontol Revy 1967;18:17-26.

31. Brannstrom M, Nyborg H. Pulp reaction to a temporary zinc oxide/eugenol cement. J Prosthet Dent 1976;35:185-91.

32. Brannstrom M, Nordenvall K-J, Torstenson B. Pulpal reaction to IRM cement: an intermediate restorative material containing eugenol. J Dent Child 1981;48(4):259-63.

33. Kanca J. An alternate hypothesis to the cause of pulpal inflammation in teeth treated with phosphoric acid on the dentin. Quintessence Int 1990;21:83-6.

34. White K, Cox C, Kanca J, Farmer J, Ramus D, Snuggs H. Histologic pulpal response of acid etching vital dentin. J Dent Res 1992;71(Special Issue):188(Abstract No. 658).

35. Willis G, Cochran M, Kafrawy A. Pulpal response to dentin acid etching in sealed and unsealed restorations. J Dent Res 1992;71(Special Issue):195(Abstract No. 717).

36. Duke E, Lindemuth J. Polymeric adhesion to dentin: contrasting substrates. Am J Dent 1990;3(6):264-70.

37. Kanca J. Dental adhesion and the All Bond system. J Esthet Dent 1991;3:129-32.

38. Souza M, Retief D, Russell C. Evaluation of the All Bond adhesive system. J Dent Res 1992;71(Special Issue):169(Abstract No. 511).



The 13th Annual

Hawaiian Dental Forum

ENJOY A VACATION...
CONTINUE YOUR EDUCATION
AN INTERNATIONAL CONTINUING EDUCATION EVENT

Maui

Jan. 30 - Feb. 6, 1993
Royal Lahaina
Resort

Kauai

February 6-13, 1993
Sheraton Kauai
Resort

- Interesting current topics presented by well credentialed speakers in morning sessions.
- Special hotel and seminar arrangements for your family, staff, and friends, including FREE buffet breakfast for resort guests.
- 18 lecture hours per week, including CPR training and certification.

The 6th Annual Humanitarian & Educational Expedition

BOLIVIA/BRAZIL

March 5-19, 1993

- Educators to provide dental training and lecture at nationally attended university programs.
- Practitioners to provide dental health care and educational services for orphans and street children in Bolivia.
- Co-hosted by local Rotary Clubs, university faculty, professional groups and the Maryknoll Missioners.
- Non-sectarian; non-political, and non-profit.



For travel information, contact:

CHIPPEWA VALLEY TRAVEL SERVICE
Chippewa Falls, Wisconsin 54729 USA
Call toll free: 1-800-826-6644 (ALL USA & CANADA)

Name _____

Address _____

City/State _____

Zip _____ Phone _____

Send information on: ☐ Hawaii ☐ Latin America.
Put me on your mailing list for next year ☐.

Sponsored by/Mail to:

Dental Seminars & Symposia, Inc.

343 Prairie View Road
Chippewa Falls, Wisconsin 54729 USA
Telephone & Fax: (715) 723-1662

VISIT US AT THE ADA ANNUAL SESSION